

- Runs through all 24 valves.
- Turns each one ON for 2 seconds before turning next valve on.
- Turns all valves off when complete or when stop is pressed.

- Must have at least one valve active

- Opens valve checkbox but only allows for one to be selected

Check All Valves

Read One Chamber

Start New Experiment

Enter New Experiment Location:

Current Valves/
Set Valves

☐ Valve 1

☐ Valve 2

☒ ...

☐ Valve 24

Select All Valves Button

Stop

Active Chamber: 1

Current CO2 :

Current H2O :

Current Temp. :

Current Flow :

Current File Location: <<Absolute Path>>

Previous Chamber: 1

Previous CO2 :

Previous H2O :

Pervious Temp. :

Previous Flow :